



Trusted for Over 70 Years

FOR IMMEDIATE RELEASE
June 26, 2019

Contact: Kelly Kramer, Marketing Director
1-800-544-3100
Kelly@cincinnati-surgical.com

Introducing the NEW Aunex Blade Remover

The Aunex blade remover makes its North American debut with Cincinnati Surgical on July 1.

This is the first blade remover to remove larger pathology blades, such as the #60 and #70 blades, in addition to standard blades size 10 to 26. Until now there has not been an option that works on larger sized blades.

The ability to remove a blade safely, without using your fingers, is key to reducing sharps injuries. The idea behind this blade remover is to easily and safely remove and store used scalpel blades for disposal. This new blade remover is ideal for use anywhere blades are removed from handles, particularly in pathology labs; dental, general practice and veterinary offices; mortuaries; universities or classroom settings.



The blade remover itself is attached to a bracket that can be wall or bench mounted. This creates a single-handed blade removal process. The design of this product also allows it to hold up to 200 blades. That's double many other blade removers on the market. Once full it can be disposed of in a proper sharps container.

This product will be available through distribution or for retail sale on our website starting July 1.

###

*For more than 70 years, Cincinnati Surgical has been an industry leader in providing cutting edge surgical products such as specialty surgical blades, handles, scalpels, suture needles, and safety products. Cincinnati Surgical offers **the sharpest, most precise, most durable blades available**. An exceptional commitment to quality products and customer service is what sets Cincinnati Surgical apart from the competition.*

Corporate Headquarters: 11331 Grooms Rd., Cincinnati Ohio 45242
Ph: 513-489-6640 | Fax: 513-489-6683 | Toll free: 1-800-544-3100 | csblades@cincinnati-surgical.com