



The specialty Swann Morton Cervical Biopsy blade is primarily used with the Long Stainless Steel Major Handle (07L) to take a “cone” biopsy from the tissue in a woman’s cervix. The tissue sample is analyzed to assess the presence of malignant and non malignant carcinoma cells.


This blade is a suitable cross for part number 379100 that has been discontinued by Beaver-Visitec. However the blade will not work with the same handle. **Only the handles listed below will fit this blade.** The cutting edge of this blade is raised to a 30 degree angle (see drawing).



Cross Reference Information

	Beaver-Visitec	Swann Morton
Cervical Biopsy Blade	379100	01C
Standard Handle	375000	07S
Long Handle	372000	07L

Technical Drawing



Swann-Morton®



Cincinnati Surgical, through our partnership with Swann Morton, brings global leading surgical blade manufacturing to the US market. From the quality of the steel to our ISO manufacturing standards and visual inspection of each blade, we are able to provide a consistent, zero defect product line that will ensure durability and consistency each time it is used. With our extensive range of products and ability to ship same day, we offer our customers a product they know they can trust when quality counts.

cincinnati-surgical.com