

Blade Remover Collection



Ensure the safe removal of surgical blades with a blade remover that meets your needs. One-handed and two-handed options are available in Cincinnati Surgical's blade remover collection.

					
Part Number	48SM*	48	48QS	48CS	48EVO
Product Name	Single-Use	Blue Box	Clicksmart BladeFlask	Aunex	Clicksmart BladeFlask EVO
Removes blades size 6 - 18 (Fitment 3)	✓	✓	✓	✓	✓
Removes blades size 20 - 27 (Fitment 4)	✓	✓	✓	✓	✓
Removes blades size 60 & 70	✗	✗	✗	✓	✓
Blade Capacity	1	100	100	50-200 varies by blade size	100
Mountable with Bracket	✗	✗	✓**	✓	✓**
Unit of Measure	Box of 100 (Case = 25)	Box of 10 (Case = 60)	Each (Case = 24)	Each (Case = 12)	Each (Case = 36)

* Available in a sterile version - part number 48SMS (50/box)

** Bracket sold separately - part number 48QSB



Raising the Standard™

Learn more at:
CincinnatiSurgical.com