



Raising the Standard™

Product Catalog



Updated July 2024



Raising the Standard™

Our Story

Since 1938, Cincinnati Surgical has been recognized as an industry leader for manufacturing high quality surgical blades, scalpels, safety products, and operating room accessories. Our manufacturing heritage in the United States (1938-1961) and Sheffield, England (1962-Present) has now come under one roof in Sheffield through a historic joint manufacturing agreement with Swann-Morton. By combining our manufacturing capabilities we are able to offer our customers the following:

- Increased range of products with continued focus on design and innovation
- In-House design and build of our own specialized machinery
- ISO 13485 Certification & In-House Sterilization
- ISO 14001 Environmental Manufacturing Standards Certification
- ISO 50001 Energy Management
- Increased capacity, production capabilities, and visual quality product inspection
- Superior customer service with same day shipping

Our Portfolio

Each brand offered has its own unique shape, size, and material composition due to historical design and manufacturing processes. There are also additional differences in the available features as laid out below and throughout this catalog.

Swann-Morton®

- 200% Visual Inspection
- Packaged with VCI Paper to protect the blade from corrosion
- Strongest/Smallest Foil Package available
- Narrow Grinding Angle resulting in increased initial keenness
- High Surface Finish
- Extended Range
- CE Mark



- 200% Visual Inspection
- Packaged with VCI Paper to protect the blade from corrosion
- Strong Foil Packaging
- Narrow Grinding Angle resulting in increased initial keenness
- Medium Surface Finish
- Standard Range

 **Lance®**

- 100% Visual Inspection
- Wide Grinding Angle resulting in high durability
- Low Surface Finish
- Limited Range
- CE Mark



Terms of Trade

To Place an Order:

Phone: 513-489-6640
Toll Free: 800-544-3100
Fax: 513-489-6683
Email: po@cincinnati-surgical.com
Mail: 11331 Grooms Rd., Cincinnati, Ohio 45242

Hours of Operation:

Monday – Thursday: 8:00am - 5:00pm EST
Friday 8:00am – 3:00pm EST

Terms:

Net 30 days with approved credit
Visa, MasterCard and American Express accepted
WE DO NOT SHIP COD

Shipping and Handling:

FOB: Cincinnati, Ohio
Unless otherwise specified, orders are shipped UPS or FedEx
Ground, prepay and add
Expedited delivery available upon request
Additional handling fees may apply

Return Policy:

Returns must have an RGA. All items must be returned in original, resalable condition within 120 days from invoice date. All returns are reviewed by our customer service department and incur a 10% restocking fee. Any non-stock items or atypical volume orders will require additional approval and may be subject to a higher restocking fee.

Distribution Policy:

A \$15 processing fee will be charged on all orders under \$100.

Minimum Order Policy:

\$25 minimum order – no distribution orders under \$25 will be accepted.

Page 1

Swann-Morton

Standard Surgical Blades
Orthopedic and General Surgery Blades
Specialty Blades
Miniature Blades
Disposable Scalpels

Page 6

Cincinnati Surgical

Standard Surgical Blades
Disposable Scalpels
Miniature Scalpels

Page 10

Lance

Standard Surgical Blades
Disposable Scalpels

Page 12

Safety

KLEEN Blade Management System by Swann-Morton
Retractable Safety Scalpels
Safety Stitch Cutter
Safety Accessories

Page 15

Stitch Cutters

Page 16

Pathology/Post Mortem

Major Surgical Blades
Disposable Knife
Post Mortem Stapler

Page 17

Skin Graft

Skin Graft Blades
Disposable Scalpel
Skin Graft Handles

Page 18

Handles

Page 21

Needles

Swann-Morton is a world leader in the manufacturing of surgical blades, scalpels and handles. It is a name respected globally for quality, precision, consistency and reliability.

Product Features:

200% Visual Inspection
Packaged with VCI Paper to protect the blade from corrosion
Strongest/Smallest Foil Package available
Narrow Grinding Angle resulting in increased initial keenness

High Surface Finish
CE Mark
Largest Range

Blades



SM Steel Blade Size 6
packaged 100/box
00SM6: sterile, carbon
01SM6: sterile, stainless
02SMV6: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size 9
packaged 100/box
00SM9: sterile, carbon
01SM9: sterile, stainless



SM Steel Blade Size 10
packaged 100/box
00SM10: sterile, carbon
01SM10: sterile, stainless
02SM10: non-sterile, carbon
02SMV10: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size 10A
packaged 100/box
00SM10A: sterile, carbon
01SM10A: sterile, stainless
02SM10A: non-sterile, carbon
02SMV10A: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size 10R
packaged 100/box
01SM10R: sterile, stainless



SM Steel Blade Size 11
packaged 100/box
00SM11: sterile, carbon
01SM11: sterile, stainless
02SM11: non-sterile, carbon
02SMV11: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size 11P
packaged 100/box
00SM11P: sterile, carbon
01SM11P: sterile, stainless
02SM11P: non-sterile, carbon
02SMV11P: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size E11
packaged 100/box
00SME11: sterile, carbon
01SME11: sterile, stainless



SM Steel Blade Size 12
packaged 100/box
00SM12: sterile, carbon
01SM12: sterile, stainless



SM Steel Blade Size 12D
packaged 100/box
00SM12D: sterile, carbon
01SM12D: sterile, stainless



SM Steel Blade Size 13
packaged 100/box
00SM13: sterile, carbon
01SM13: sterile, stainless



SM Steel Blade Size 14
packaged 100/box
00SM14: sterile, carbon
01SM14: sterile, stainless
02SMV15C: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size 15
packaged 100/box
00SM15: sterile, carbon
01SM15: sterile, stainless
02SM15: non-sterile, carbon
02SMV15: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size 15A
packaged 100/box
00SM15A: sterile, carbon
01SM15A: sterile, stainless



SM Steel Blade Size 15C
packaged 100/box
00SM15C: sterile, carbon
01SM15C: sterile, stainless
02SMV15C: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size 15T
packaged 100/box
00SM15T: sterile, carbon
01SM15T: sterile, stainless

*Blades are not shown to actual size.



SM Steel Blade Size 16
packaged 100/box
00SM16: sterile, carbon
01SM16: sterile, stainless
02SMV16: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size 18
packaged 100/box
00SM18: sterile, carbon
01SM18: sterile, stainless
02SMV18: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size 19
packaged 100/box
01SM19: sterile, stainless



SM Steel Blade Size 20
packaged 100/box
00SM20: sterile, carbon
01SM20: sterile, stainless
02SMV20: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size 21
packaged 100/box
00SM21: sterile, carbon
01SM21: sterile, stainless
02SM21: non-sterile, carbon



SM Steel Blade Size 22
packaged 100/box
00SM22: sterile, carbon
01SM22: sterile, stainless
02SM22: non-sterile, carbon



SM Steel Blade Size 22A
packaged 100/box
00SM22A: sterile, carbon
01SM22A: sterile, stainless



SM Steel Blade Size 23
packaged 100/box
00SM23: sterile, carbon
01SM23: sterile, stainless
02SM23: non-sterile, carbon
02SMV23: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size 24
packaged 100/box
00SM24: sterile, carbon
01SM24: sterile, stainless
02SM24: non-sterile, carbon
02SMV24: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size 25
packaged 100/box
00SM25: sterile, carbon
01SM25: sterile, stainless
02SMV25: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size 25A
packaged 100/box
00SM25A: sterile, carbon
01SM25A: sterile, stainless
02SMV25A: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size 26
packaged 100/box
00SM26: sterile, carbon
01SM26: sterile, stainless
02SMV26: non-sterile carbon
(20 PK/5 VPI)



SM Steel Blade Size 27
packaged 100/box
00SM27: sterile, carbon
01SM27: sterile, stainless



SM Steel Blade Size 36
packaged 100/box
00SM36: sterile, carbon



SM Steel Blade Size 40
packaged 100/box
00SM40: sterile, carbon



Cervical Biopsy Blade
packaged 10/box
01C: sterile, stainless



Myringotomy Blade
packaged 10/box
01M: sterile, stainless



Podiatry Blade
packaged 5/box
PD81: non-sterile, stainless



Podiatry Blade
packaged 5/box
PD82: non-sterile, stainless



Miniature Surgical Blade Size 61
packaged 25/box
01SM61: sterile, stainless



Miniature Surgical Blade Size 62
packaged 25/box
01SM62: sterile, stainless



Miniature Surgical Blade Size 63
packaged 25/box
01SM63: sterile, stainless



Miniature Surgical Blade Size 64
packaged 25/box
01SM64: sterile, stainless



Miniature Surgical Blade Size 65
packaged 25/box
01SM65: sterile, stainless



Miniature Surgical Blade Size 65A
packaged 25/box
01SM65A: sterile, stainless



Miniature Surgical Blade Size 67
packaged 25/box
01SM67: sterile, stainless



Miniature Surgical Blade Size 68
packaged 25/box
01SM68: sterile, stainless



Miniature Surgical Blade Size 69
packaged 25/box
01SM69: sterile, stainless



Miniature Surgical Blade Size 90
packaged 25/box
01SP90: sterile, stainless




Miniature Surgical Blade Size 91
packaged 25/box
01SP91: sterile, stainless

*Blades are not shown to actual size.

*Scalpels are not shown to actual size.

SM Disposable Scalpel Size 9


SM0517: sterile, stainless
packaged 10/box



SM Disposable Scalpel Size 10

SM0501: sterile, stainless
packaged 10/box

SM2331: non-sterile, stainless
packaged 1000/box




SM Disposable Scalpel Size 10R

SM6690: sterile, stainless
packaged 10/box



SM Disposable Scalpel Size 11

SM0503: sterile, stainless
packaged 10/box



SM Disposable Scalpel Size 11P


SM0591: sterile, stainless
packaged 10/box

SM2391: non-sterile, stainless
packaged 1000/box

SM Disposable Scalpel Size 12


SM0504: sterile, stainless
packaged 10/box

SM2334: non-sterile, stainless
packaged 1000/box




SM Disposable Scalpel Size 12D

SM0518: sterile, stainless
packaged 10/box



SM Disposable Scalpel Size 14

SM0519: sterile, stainless
packaged 10/box




*Scalpels are not shown to actual size.

SM Disposable Scalpel Size 15


SM0505: sterile, stainless
packaged 10/box

SM2335: non-sterile, stainless
packaged 1000/box




SM Disposable Scalpel Size 15A

SM0520: sterile, stainless
packaged 10/box



SM Disposable Scalpel Size 15T


SM0592: sterile, stainless
packaged 10/box



SM Disposable Scalpel Size 20

SM0506: sterile, stainless
packaged 10/box


SM2336: non-sterile, stainless
packaged 1000/box



SM Disposable Scalpel Size 21


SM0507: sterile, stainless
packaged 10/box

SM2341: non-sterile, stainless
packaged 1000/box




SM Disposable Scalpel Size 22

SM0508: sterile, stainless
packaged 10/box




SM Disposable Scalpel Size 23

SM0510: sterile, stainless
packaged 10/box



SM Disposable Double Edged SGD Scalpel

SMSGD: sterile, stainless
packaged 10/box



Disposable Mini Scalpel Size 10

packaged 100/box

1110SM: Swann Morton
non-sterile, stainless



Disposable Mini Scalpel Size 11P

packaged 100/box

1111SM: Swann Morton
non-sterile, stainless



Disposable Mini Scalpel Size 15

packaged 100/box

1115SM: Swann Morton
non-sterile, stainless



*Scalpels are not shown to actual size.



Raising the Standard™

Cincinnati Surgical is a leading manufacturer of high grade British surgical steel products. For more than 70 years we have provided our customers with quality products and a high level of customer service. Cincinnati Surgical is also known as one of the original manufacturers of stainless steel blades.

Product Features:

200% Visual Inspection

Packaged with VCI Paper to protect the blade from corrosion

Narrow Grinding Angle resulting in increased initial keenness

Medium Surface Finish

Strong Foil Package

Blades



CS Steel Blade Size 10
packaged 100/box
0010: sterile, carbon
0110: sterile, stainless
0210: non-sterile, carbon
0410: non-sterile stainless
(500/pkg)



CS Steel Blade Size 11
packaged 100/box
0011: sterile, carbon
0111: sterile, stainless
0211: non-sterile, carbon
0411: non-sterile stainless
(500/pkg)



CS Steel Blade Size 11P
packaged 500/pkg
0411P: non-sterile, stainless



CS Steel Blade Size 12
packaged 100/box
0012: sterile, carbon
0112: sterile, stainless
0212: non-sterile, carbon
0412: non-sterile stainless
(500/pkg)



CS Steel Blade Size 12D
packaged 100/box
0112D: sterile, stainless



CS Steel Blade Size 15
packaged 100/box
0015: sterile, carbon
0115: sterile, stainless
0215: non-sterile, carbon
0415: non-sterile stainless
(500/pkg)



CS Steel Blade Size 15C
packaged 100/box
0115C: sterile, stainless



CS Steel Blade Size 20
packaged 100/box
0020: sterile, carbon
0120: sterile, stainless
0220: non-sterile, carbon



CS Steel Blade Size 21
packaged 100/box
0021: sterile, carbon
0121: sterile, stainless
0221: non-sterile, carbon
0421: non-sterile stainless
(500/pkg)



CS Steel Blade Size 22
packaged 100/box
0022: sterile, carbon
0122: sterile, stainless
0222: non-sterile, carbon
0422: non-sterile stainless
(500/pkg)



CS Steel Blade Size 23
packaged 100/box
0023: sterile, stainless
0123: sterile, stainless
0223: non-sterile, carbon
0423: non-sterile stainless
(500/pkg)



CS Steel Blade Size 24
packaged 100/box
0224: non-sterile, carbon



CS Steel Blade Size 25A
packaged 100/box
0225A: non-sterile carbon

*Blades are not shown to actual size.

CS Disposable Scalpel Size 10

0510: sterile, stainless 10/box

0610: non-sterile, stainless
100/box



CS Disposable Scalpel Size 11

0511: sterile, stainless 10/box

0611: non-sterile, stainless
100/box



CS Disposable Scalpel Size 12

0512: sterile, stainless 10/box

0612: non-sterile, stainless
100/box



CS Disposable Scalpel Size 15

0515: sterile, stainless 10/box

0615: non-sterile, stainless
100/box



CS Disposable Scalpel Size 20

0520: sterile, stainless 10/box

0620: non-sterile, stainless
100/box



CS Disposable Scalpel Size 21

0521: sterile, stainless 10/box

0621: non-sterile, stainless
100/box



CS Disposable Scalpel Size 22

0522: sterile, stainless 10/box

0622: non-sterile, stainless
100/box



*Scalpels are not shown to actual size.

Disposable Mini Scalpel Size 10
packaged 100/box

1110SM:Swann Morton
non-sterile, stainless



Disposable Mini Scalpel Size 11P
packaged 100/box

1111:Cincinnati Surgical
non-sterile, stainless



Disposable Mini Scalpel Size 15
packaged 100/box

1115SM:Swann Morton
non-sterile, stainless





Lance products enjoy a reputation of long lasting cutting edges. Processed from British “cold rolled” steel which is superior in tolerance, concentricity, and straightness allowing for wider grinding angles.

Product Features:

100% Visual Inspection
Wide Grinding Angle resulting in high durability
Low Surface Finish

Limited Range
CE Mark

Blades



Lance Steel Blade Size 10
packaged 100/box
92010: sterile, carbon
92110: sterile, stainless



Lance Steel Blade Size 11
packaged 100/box
92011: sterile, carbon
92111: sterile, stainless



Lance Steel Blade Size 12
packaged 100/box
92012: sterile, carbon
92112: sterile, stainless



Lance Steel Blade Size 15
packaged 100/box
92015: sterile, carbon
92115: sterile, stainless



Lance Steel Blade Size 20
packaged 100/box
92020: sterile, carbon
92120: sterile, stainless



Lance Steel Blade Size 21
packaged 100/box
92021: sterile, carbon
92121: sterile, stainless



Lance Steel Blade Size 22
packaged 100/box
92022: sterile, carbon
92122: sterile, stainless



Lance Steel Blade Size 23
packaged 100/box
92023: sterile, carbon

Lance Disposable Scalpel Size 10

92410: sterile, stainless 10/box



Lance Disposable Scalpel Size 11


92411: sterile, stainless 10/box

92711: non-sterile, stainless
100/box



Lance Disposable Scalpel Size 12

92412: sterile, stainless 10/box



Lance Disposable Scalpel Size 15

92415: sterile, stainless 10/box

92715: non-sterile, stainless
100/box



Lance Disposable Scalpel Size 20

92420: sterile, stainless 10/box




Lance Disposable Scalpel Size 21

92421: sterile, stainless 10/box



Lance Disposable Scalpel Size 22

92422: sterile, stainless 10/box



Lance Disposable Scalpel Size 23

92423: sterile, stainless 10/box



KLEEN Blade Management System by *Swann-Morton®*




SM Steel Blade Size 10

01SM10K:

sterile, stainless 50/box

03SM10K:

non-sterile, stainless 500/box




SM Steel Blade Size 11

01SM11K:

sterile, stainless 50/box

03SM11K:

non-sterile, stainless 500/box




SM Steel Blade Size 11P

01SM11PK:

sterile, stainless 50/box

03SM11PK:

non-sterile, stainless 500/box




SM Steel Blade Size 15

01SM15K:

sterile, stainless 50/box

03SM15K:

non-sterile, stainless 500/box




SM Steel Blade Size 20

01SM20K:

sterile, stainless 50/box

03SM20K:

non-sterile, stainless 500/box



SM Steel Blade Size 21

01SM21K:

sterile, stainless 50/box

03SM21K:

non-sterile, stainless 500/box




SM Steel Blade Size 22

01SM22K:

sterile, stainless 50/box

03SM22K:

non-sterile, stainless 500/box



SM Steel Blade Size 23

01SM23K:

sterile, stainless 50/box

03SM23K:

non-sterile, stainless 500/box



SM Steel Blade Size 24

01SM24K:

sterile, stainless 50/box

03SM24K:

non-sterile, stainless 500/box

Safety Scalpels

CS Safety Scalpel Size 10

05S10: sterile, stainless 25/box

06S10: non-sterile, stainless 50/box

SM Safety Scalpel Size 10

05SM10: sterile, stainless 25/box

06SM10: non-sterile, stainless 500/box



CS Safety Scalpel Size 11P

05S11: sterile, stainless 25/box

06S11: non-sterile, stainless 50/box

SM Safety Scalpel Size 11P

05SM11P: sterile, stainless 25/box

06SM11P: non-sterile, stainless 500/box



CS Safety Scalpel Size 15

05S15: sterile, stainless 25/box

06S15: non-sterile, stainless 50/box

SM Safety Scalpel Size 15

05SM15: sterile, stainless 25/box

06SM15: non-sterile, stainless 500/box



*Items are not shown to actual size.

SM Safety Scalpel Size 20**05SM20:** sterile, stainless 25/box**06SM20:** non-sterile, stainless 500/box**SM Safety Scalpel Size 21****05SM21:** sterile, stainless 25/box**06SM21:** non-sterile, stainless 500/box**SM Safety Scalpel Size 22****05SM22:** sterile, stainless 25/box**06SM22:** non-sterile, stainless 500/box**SM Safety Scalpel Size 23****05SM23:** sterile, stainless 25/box**06SM23:** non-sterile, stainless 500/box**SM Safety Scalpel Size 24****05SM24:** sterile, stainless 25/box**06SM24:** non-sterile, stainless 500/box**CS Safety Scalpel Size 25A****05S25:** sterile, stainless 25/box**06S25:** non-sterile, stainless 500/box

Safety Stitch Cutter

Safety Stitch Cutter**SM3926:** sterile, stainless 25/box**SM4026:** non-sterile, stainless 500/box

*Items are not shown to actual size.



Clicksmart BladeFlask
48QS: packaged individually

Removes standard sized surgical blades size 6 - 27



Clicksmart BladeFlask EVO
48EVO: packaged individually

Removes all standard sized surgical blades including size 60 and 70 blades



Clicksmart BladeFlask Bracket
48QSB: packaged individually

For Use on both version of the BladeFlask blade remover.

Aunex Blade Remover
48CS: packaged individually

Removes all standard sized surgical blades including size 60 and 70 blades

Holds up to 200 blades (This number will vary based on the size of the blade being removed.)




Cincinnati Surgical Blade Remover
48: packaged 10/box



Fits 100 blades


Single Use Blade Remover
48SM: Non-sterile, packaged 100/box
48SMS: Sterile, packaged 50/box



*Items are not shown to actual size.

SM Disposable Stitch Cutter Scalpel

SM0526: sterile, stainless 10/box



CS Stitch Cutter

44: sterile, stainless 100/box



SM Stitch Cutter Blade

SM0326: sterile, stainless blade 100/box
Fits any #3 fitment handle




SM Stitch Cutter

SM0420: sterile, carbon, short 100/box



SM Stitch Cutter

SM0421: sterile, stainless, long 100/box



SM Stitch Cutter

SM0422: sterile, stainless, midi 100/box



Safety Stitch Cutter

SM3926: sterile, stainless 25/box

SM4026: non-sterile, stainless 500/box



*Items are not shown to actual size.



SM Surgical Blade Size 22B
packaged 100/box
02SM22B: non-sterile, carbon



CS Surgical Blade Size 60
packaged 100/box
0260: non-sterile, carbon
0260CS: non-sterile, carbon
(single foil)
Fits Handles: 0740S, 0780S
0780SM-S



SM Surgical Blade Size 60
packaged 50/box
02SM60: non-sterile, carbon
Fits Handles: 0740S, 0780S
0780SM-S



SM Surgical Blade Size 60B
packaged 50/box
02SM60B: non-sterile, carbon
Fits Handles: 0740S, 0780S
0780SM-S



CS Surgical Blade Size 70
packaged 100/box
0270: non-sterile, carbon

Fits Handles: 0740S, 0780S
0780SM



SM Surgical Blade Size 40
packaged 10/box
02PM40: non-sterile, carbon

Fits Handles: 07PM, 07PMH



SM Surgical Blade Size 40B
packaged 10/box
02PM40B: non-sterile, carbon

Fits Handles: 07PM, 07PMH

Disposable Knife Size 40
packaged individually
06PM40: non-sterile, carbon



Post Mortem Stapler
packaged individually or 50/case
853590: non-sterile
pre-loaded with 35 staples



Features include an easy to read staple counter, an arrow guide for accurate placement and a low profile for a clear view of the work area. It is not for medical use.

*Items are not shown to actual size.

Skin Graft Blade SG3

packaged 20/box

SM2201: sterile, stainless

Fits Handles:

073S, 073L



Skin Graft Blade SGD

packaged 20/box

SM2203: sterile, stainless

Fits Handles:

073S, 073L



Skin Graft Blade 158MM

packaged 10/box

SM9940: sterile, stainless

Fits Handles:

SM9901, SM9902,

SM9903



Skin Graft Blade Silvers

packaged 10/box

SM9942: sterile, stainless

Fits Handle:


SM9911



SM Disposable Double Edged SGD Scalpel

packaged 10/box

SM9940: sterile, stainless



SM Braithwaite Handle

SM9901: stainless handle

SM9904: stainless handle with
2 boxes of SM9940



SM Cobbett Handle


SM9902: stainless handle



SM Watson Handle

SM9903: stainless handle

SM9906: stainless handle with
2 boxes of SM9940



SM Silvers Handle

SM9911: stainless handle



*Items are not shown to actual size.

#3 Fitment Handles

Handle Size 3

07B3: stainless, lightweight
Fits Blades 6 - 16



Handle Size 3

073S: stainless
Fits Blades 6 - 16
UDI Compliant



Handle Size 3

073S-NU: stainless
Fits Blades 6 - 16



Handle Size 3 Long

073L: stainless
Fits Blades 6 - 16



Handle Size 3 Slim

07SM9: stainless
Fits Blades 6 - 16



Handle Size 5

07SM5B-S: stainless
Fits Blades: 6 - 16



Handle Size 5

07SM5B-E: stainless
Fits Blades: 6 - 16



Handle Size 7

077S: stainless
Fits Blades: 6 - 16



*Handles are not shown to actual size.

#4 Fitment Handles

Handle Size 4

074S: stainless
Fits Blades 20-26
UDI Compliant



Handle Size 4

074S-NU: stainless
Fits Blades 20-26



Handle Size 4 Long

074L: stainless
Fits Blades 20-26



Handle Size 6

07SM6B-S: stainless
Fits Blades: 20 - 26



Handle Size 6

07SM6B-E: ergonomic
Fits Blades: 20 - 26



Handle Size 8

078S: stainless
Fits Blades: 20 - 26



Major Handle

07L: stainless, long
Fits Blade: 01C



Major Handle

07S: stainless, standard
Fits Blade: 01C



*Handles are not shown to actual size.

Miniature Handle Size 13cm**07SF13:** stainless, quick release feature

Fits Blades: 01M and 61-91

**Miniature Handle Size 10cm****07SF23:** stainless, quick release feature

Fits Blades: 01M and 61-91

**Miniature Handle Size 7.5cm****07SF3:** stainless

Fits Blades: 01M and 61-91

**Miniature Handle Size 5cm****07SF4:** stainless

Fits Blades: 01M and 61-91

**Handle****0780SM-E:** ergonomic, pathology/post mortem


Fits Blades: 0260, 0260CS, 02SM60

02SM60B and 0270

**Handle****0780SM-S:** stainless, pathology/post mortem

Fits Blades: 0260, 0260CS, 02SM60

02SM60B and 0270

**Handle Size PM40****07PM:** stainless, pathology/post mortem

Fits Blades: PM40 and PM40B

**Handle Size PM40****07PMH:** stainless, pathology/post mortem

Fits Blades: PM40 and PM40B



*Handles are not shown to actual size.

Regular Surgeon Needles

3/8 Circle Cutting Edge

Spring Eye

08622: non-sterile, stainless

12/pkg

Available Sizes	Dimensions
12	.69mm x 30mm
18	.61mm x 22mm




Post Mortem Needles

3/8 Circle Cutting Edge

08816: non-sterile, stainless

12/pkg

Available Sizes	Dimensions
3	1.42mm x 75mm
4	1.80mm x 100mm
5	2.00mm x 125mm
6	2.00mm x 152mm




Straight Cutting Needles

Triangular Point

08817: non-sterile, stainless

12/pkg

Available Sizes	Dimensions
2	1.02mm x 2"
3	1.27mm x 3"
4	1.42mm x 4"




Post Mortem Needles

Serpentine Double Curved Cutting Edge

08819: non-sterile, stainless

12/pkg

Available Sizes	Dimensions
0	5"
1	4.75"
2	4.25"
3	3.25"
4	2.25"




Post Mortem Needles

Half Curved Cutting Edge

08820: non-sterile, stainless

12/pkg

Available Sizes	Dimensions
00	2.20mm x 150mm
0	2.20mm x 140mm
1	2.00mm x 114mm
2	1.8mm x 89mm




*Needles are not shown to actual size.

Regular Surgeon Needles Half Curved Cutting Edge

08821: non-sterile, stainless
12/pkg


Available Sizes	Dimensions
1	1.21mm x 80mm
2	1.20mm x 68mm
4	1.20mm x 64mm
6	1.10mm x 55mm
8	1.00mm x 50mm
10	.93mm x 43mm
12	.76mm x 35mm
14	.69mm x 29mm
16	.60mm x 27mm
18	.64mm x 25mm
20	.53mm x 22mm
22	.50mm x 19mm



Regular Surgeon Needles 3/8 Circle Cutting Edge

08822: non-sterile, stainless
12/pkg



Available Sizes	Dimensions
1	1.20mm x 80mm
2	1.22mm x 68mm
4	1.21mm x 64mm
6	1.10mm x 55mm
8	1.00mm x 50mm
10	.94mm x 43mm
12	.76mm x 35mm
14	.69mm x 29mm
16	.64mm x 27mm
18	.53mm x 22mm
20	.48mm x 17mm



Regular Surgeon Needles 3/8 Circle Taper Point

08822T: non-sterile, stainless
12/pkg

Available Sizes	Dimensions
4	1.16mm x 65mm
6	1.09mm x 55mm
12	.76mm x 35mm
14	.69mm x 29mm
16	.64mm x 27mm
18	.53mm x 22mm




Mayo's Catgut Needles 1/2 Circle Taper Point

08824: non-sterile, stainless
12/pkg

Available Sizes	Dimensions
2	1.40mm x 45mm
3	1.40mm x 38mm
4	1.40mm x 35mm
6	1.30mm x 28mm

09824: sterile, stainless
50/box


Available Sizes	Dimensions
2	1.40mm x 45mm
3	1.40mm x 38mm
4	1.40mm x 35mm
5	1.30mm x 32mm
6	1.30mm x 28mm



Mayo's Catgut Needles 1/2 Circle Trocar Point

09826: sterile, stainless
50/box

Available Sizes	Dimensions
2	1.40mm x 45mm
3	1.40mm x 38mm
4	1.40mm x 35mm
5	1.30mm x 32mm
6	1.30mm x 28mm



*Needles are not shown to actual size.


Keith's Straight Abdominal Needles Triangular Point

08827: non-sterile, stainless
12/pkg

Available Sizes	Dimensions
1 ½	.74mm x 38mm
1 ¾	.75mm x 44mm
2	.76mm x 51mm
2 ½	.81mm x 64mm
2 7/8	.92mm x 76mm
3 7/8	1.00mm x 95mm

09827: sterile, stainless
50/box

Available Sizes	Dimensions
2	.76mm x 51mm
2 ½	.81mm x 64mm
2 7/8	.83mm x 76mm



Eye Needles 3/8 Circle Cutting Edge

08830: non-sterile, stainless
12/pkg


Available Sizes	Dimensions
1	.53mm x 28mm
2	.53mm x 22mm



Eye Needles 1/2 Circle Cutting Edge

08831: non-sterile, stainless
12/pkg


Available Sizes	Dimensions
1	.53mm x 25mm
2	.53mm x 22mm



Fistula Needles 1/2 Circle Cutting Edge

08832: non-sterile, stainless
12/pkg


Available Sizes	Dimensions
0	1.42mm x 64mm
9	.76mm x 19mm



Regular Surgeon Needles 1/2 Circle Cutting Edge

08834: non-sterile, stainless
12/pkg

Available Sizes	Dimensions
00	1.20mm x 68mm
0	1.20mm x 64mm
1	1.09mm x 59mm
2	1.00mm x 50mm
4	.90mm x 43mm
6	.68mm x 29mm
9	.50mm x 19mm




*Needles are not shown to actual size.

Regular Surgeon Needles 1/2 Circle Taper Point

08834T: non-sterile, stainless
12/pkg


Available Sizes	Dimensions
0	1.20mm x 64mm
2	1.00mm x 50mm
4	.90mm x 43mm
6	.68mm x 29mm
9	.53mm x 19mm



Mayo's Intestinal Needles 3/8 Circle Taper Point

08837: non-sterile, stainless
12/pkg

Available Sizes	Dimensions
1	.53mm x 25mm
1 3/4	.76mm x 45mm




Ferguson's Needles 1/2 Circle Taper Point

08842: non-sterile, stainless
12/pkg

Available Sizes	Dimensions
4	1.22mm x 60mm
6	1.10mm x 55mm
8	.92mm x 44mm
10	.76mm x 33mm
18	.53mm x 19mm
20	.64mm x 15mm
22	.53mm x 13mm

09842: sterile, stainless
50/box

Available Sizes	Dimensions
2	1.22 mm x 69mm



Lane's Cleft Palate Needles 1/2 Circle Cutting Edge

08844: non-sterile, stainless
12/pkg


Available Sizes	Dimensions
1	.71mm x 25mm



Loopuyt's Cutting Edge Needles

08846: non-sterile, stainless
12/pkg

Available Sizes	Dimensions
1	1.42mm x 80mm
2	1.22mm x 68mm
3	1.10mm x 60mm
4	1.00mm x 52mm
5	1.00mm x 46mm



*Needles are not shown to actual size.

Lane's Cleft Palate Needles 3/8 Circle Cutting Edge

08848: non-sterile, stainless
12/pkg

Available Sizes	Dimensions
1	.71mm x 22mm
2	.69mm x 19mm



Davis Tonsil Needles 1/2 Circle Taper Point

09849: sterile, stainless
50/box


Available Sizes	Dimensions
one size	1.20mm x 25mm



Martin's Uterine Needles 1/2 Circle Cutting Edge

08860: non-sterile, stainless
12/pkg


Available Sizes	Dimensions
1	1.60mm x 74mm
3	1.27mm x 54mm
8	.86mm x 22mm



French Spring Eye Needles 1/2 Circle Taper Point

09871: sterile, stainless
50/box


Available Sizes	Dimensions
2	.76mm x 26mm



Regular Surgeon Needles 3/8 Circle Reverse Cutting Edge

08922: non-sterile, stainless
12/pkg


Available Sizes	Dimensions
6	1.10mm x 55mm
8	1.00mm x 50mm
10	.93mm x 43mm
12	.76mm x 35mm
14	.69mm x 29mm
16	.64mm x 27mm
18	.53mm x 22mm



Regular Surgeon Needles 1/2 Circle Reverse Cutting Edge

08934: non-sterile, stainless
12/pkg

Available Sizes	Dimensions
3	1.00mm x 50mm
5	.85mm x 38mm
7	.60mm x 27mm





Raising the Standard™

THE SHARPEST, MOST PRECISE,
MOST DURABLE BLADES FOR OVER

80 YEARS

CincinnatiSurgical.com